


# LONRF1 Antibody / RNF191 (RQ8547)

Catalog No.	Formulation	Size
RQ8547	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

## **Bulk quote request**

Availability	1-3 days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit IgG
Purity	Antigen affinity purified
Buffer	Lyophilized from 1X PBS with 2% Trehalose
UniProt	Q17RB8
Applications	Western Blot : 0.5-1ug/ml ELISA : 0.1-0.5ug/ml
Limitations	This LONRF1 antibody is available for research use only.



Western blot testing of human 1) HeLa, 2) RT4, 3) HEL and 4) HepG2 cell lysate with LONRF1 antibody. Predicted molecular weight: 85-87 kDa.

## **Description**

LONRF1, also known as LON peptidase N-terminal domain and RING finger protein 1 (RNF191), consists of an 86.7 kDa and an 85.5 kDa isoform. There is no information regarding the function of the protein or the research being conducted with the protein. The protein interacts with TRIM9, UBE2L6, and BRCA1 proteins.

#### **Application Notes**

Optimal dilution of the LONRF1 antibody should be determined by the researcher.

## **Immunogen**

An E.coli-derived human recombinant protein (amino acids L52-D442) was used as the immunogen for the LONRF1 antibody.

# **Storage**

After reconstitution, the LONRF1 antibody can be stored for up to one month at 4oC. For long-term, aliquot and store at -20oC. Avoid repeated freezing and thawing.