


LRP6 Antibody (R33281)

Catalog No.	Formulation	Size
R33281-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	4040
Applications	Immunohistochemistry (FFPE) : 4-6ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This LRP6 antibody is available for research use only.



IHC testing of FFPE human kidney with LRP6 antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, HRP-staining.

Description

Additional name(s) for this target protein: Low density lipoprotein receptor-related protein 6

Application Notes

Optimal dilution of the LRP6 antibody should be determined by the researcher.

Immunogen

Amino acids CMTSVATAKGYTSDL were used as the immunogen for this LRP6 antibody.

Storage

Aliquot and store the LRP6 antibody at -20oC.