


LRRK1 Antibody (F40142)

Catalog No.	Formulation	Size
F40142-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F40142-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Format	Purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Purified
UniProt	Q38SD2
Applications	Western Blot : 1:1000
Limitations	This LRRK1 antibody is available for research use only.



Description

Protein kinases are enzymes that transfer a phosphate group from a phosphate donor, generally the γ phosphate of ATP, onto an acceptor amino acid in a substrate protein. By this basic mechanism, protein kinases mediate most of the signal transduction in eukaryotic cells, regulating cellular metabolism, transcription, cell cycle progression, cytoskeletal rearrangement and cell movement, apoptosis, and differentiation. With more than 500 gene products, the protein kinase family is one of the largest families of proteins in eukaryotes. The family has been classified in 8 major groups based on sequence comparison of their tyrosine (PTK) or serine/threonine (STK) kinase catalytic domains.

Application Notes

Titration of the LRRK1 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 1981-2015 from the human protein was used as the immunogen for this LRRK1 antibody.

Storage

Aliquot the LRRK1 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.