


M6prbp1 Antibody (R34889)

Catalog No.	Formulation	Size
R34889-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	66905
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This M6prbp1 antibody is available for research use only.



Western blot testing of rat skin lysate with M6prbp1 antibody at 0.1ug/ml. Predicted molecular weight: ~47kDa.

Description

Additional name(s) for this target protein: Mannose-6-phosphate receptor binding protein 1; Plin3; Perilipin 3; TIP47

Application Notes

Optimal dilution of the M6prbp1 antibody should be determined by the researcher.

Immunogen

Amino acids NQKTPQDAEKDPAK were used as the immunogen for this M6prbp1 antibody.

Storage

Aliquot and store the M6prbp1 antibody at -20oC.