


## MAGOH Antibody (R35962)

Catalog No.	Formulation	Size
R35962-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	4116
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:2000
Limitations	This MAGOH antibody is available for research use only.



Western blot testing of nuclear HeLa lysate with MAGOH antibody at 1ug/ml. Predicted molecular weight: ~17kDa.

### Description

Additional name(s) for this target protein: Mago-nashi homolog, MAGOHA, MAGOH

### Application Notes

Optimal dilution of the MAGOH antibody should be determined by the researcher.

## Immunogen

Amino acids SLIGLHFKIKPI were used as the immunogen for this MAGOH antibody.

## Storage

Aliquot and store the MAGOH antibody at -20oC.