


Malt1 Antibody (isoform 2) (R35705)

Catalog No.	Formulation	Size
R35705-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	240354
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This Malt1 antibody is available for research use only.



Western blot testing of mouse kidney lysate with Malt1 antibody at 2ug/ml. Predicted molecular weight: ~92kDa.

Description

Additional name(s) for this target protein: Mucosa associated lymphoid tissue lymphoma translocation gene 1

Application Notes

Optimal dilution of the Malt1 antibody should be determined by the researcher.

1. This Malt1 antibody is specific to isoform 2.

Immunogen

Amino acids RPDNKEQTGQPLAK were used as the immunogen for this Malt1 antibody.

Storage

Aliquot and store the Malt1 antibody at -20oC.