


## MARCH10 Antibody (R36050)


Catalog No.	Formulation	Size
R36050-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Rat
Predicted Reactivity	Cow
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	162333
Applications	Western Blot : 0.3-1ug/ml IHC (FFPE) : 5ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This MARCH10 antibody is available for research use only.



Western blot testing of rat testis lysate with MARCH10 antibody at 0.3ug/ml. Predicted/observed molecular weight: ~90 kDa. Both observed bands are blocked by addition of immunizing peptide.



IHC testing of FFPE human testis tissue with MARCH10 antibody at 5ug/ml. Required HIER: steamed antigen retrieval with pH6 citrate buffer; AP-staining.

## Description

Additional name(s) for this target protein: Membrane-associated ring finger (C3HC4) 10, MARCH-X

## Application Notes

Optimal dilution of the MARCH10 antibody should be determined by the researcher.

## Immunogen

Amino acids SRNYPQPRTEE were used as the immunogen for this MARCH10 antibody.

## Storage

Aliquot and store the MARCH10 antibody at -20oC.