


# MARCH6 Antibody (R35284)

Catalog No.	Formulation	Size
R35284-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide $$	100 ug

#### **Bulk quote request**

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	10299
Applications	Western Blot: 0.2-0.5ug/ml ELISA (peptide) LOD: 1:16000
Limitations	This MARCH6 antibody is available for research use only.



### **Description**

Additional name(s) for this target protein: Membrane-associated ring finger (C3HC4) 6, TEB4, MARCH-VI

## **Application Notes**

Optimal dilution of the MARCH6 antibody should be determined by the researcher.

#### **Immunogen**

Amino acids QGRRVIFQKVKE were used as the immunogen for this MARCH6 antibody. **Storage** Aliquot and store the MARCH6 antibody at -20oC.