


MARK4 Antibody (R34335)

Catalog No.	Formulation	Size
R34335-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	57787
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This MARK4 antibody is available for research use only.



Western blot testing of human prostate lysate with MARK4 antibody at 0.5ug/ml.
Predicted molecular weight: ~75kDa.

Description

Additional name(s) for this target protein: MAP/microtubule affinity-regulating kinase 4

Application Notes

Optimal dilution of the MARK4 antibody should be determined by the researcher.

1. This MARK4 antibody is specific for isoform 1 (long form).

Immunogen

Amino acids IKRKPPPGCSDSPGV were used as the immunogen for this MARK4 antibody.

Storage

Aliquot and store the MARK4 antibody at -20oC.