


MAS1L Antibody / MAS1-like / Mas-related G-protein coupled receptor MRG [clone 29M28] (RQ8774)

Catalog No.	Formulation	Size
RQ8774	Antibody in PBS with 0.02% sodium azide, 50% glycerol and 0.4-0.5mg/ml BSA	100 ul


Recombinant **RABBIT MONOCLONAL**

[Bulk quote request](#)

Availability	1-3 days
Species Reactivity	Human
Format	Purified
Host	Rabbit
Clonality	Recombinant Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	29M28
Purity	Affinity chromatography
UniProt	P35410
Applications	Western Blot : 1:500-1:2000 Immunoprecipitation : 1:50
Limitations	This MAS1L antibody is available for research use only.



Immunoprecipitation of MAS1L protein from 500ug of human HeLa whole cell lysate with MAS1L antibody at 1:50. Western blot testing done at 1:1000.



Western blot testing of human HeLa cell lysate with MAS1L antibody. Predicted molecular weight ~42 kDa.

Description

Mas-related G-protein coupled receptor MRG is a protein that in humans is encoded by the MAS1L gene.

Application Notes

Optimal dilution of the MAS1L antibody should be determined by the researcher.

Immunogen

A synthetic peptide specific to human MAS1L protein was used as the immunogen for the MAS1L Antibody.

Storage

Store the MAS1L antibody at -20oC.