


## MASK Antibody / STK26 / MST4 (R34517)


Catalog No.	Formulation	Size
R34517-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

<b>Availability</b>	1-3 business days
<b>Species Reactivity</b>	Human
<b>Predicted Reactivity</b>	Mouse, Rat, Dog, Pig, Cow
<b>Format</b>	Antigen affinity purified
<b>Host</b>	Goat
<b>Clonality</b>	Polyclonal (goat origin)
<b>Isotype</b>	Goat Ig
<b>Purity</b>	Antigen affinity
<b>UniProt</b>	Q9P289
<b>Gene ID</b>	51765
<b>Localization</b>	Cytoplasmic
<b>Applications</b>	Western Blot : 0.5-1.5ug/ml IHC (FFPE) : 2.5ug/ml ELISA (peptide) LOD : 1:32000
<b>Limitations</b>	This MASK antibody is available for research use only.



Western blot testing of Jurkat lysate with MASK antibody at 0.5ug/ml. Predicted molecular weight: ~47 kDa.



IHC testing of FFPE human thymus tissue with MASK antibody at 2.5ug/ml. Required HIER: steamed antigen retrieval with pH6 citrate buffer; AP-staining.

## Description

Additional name(s) for this target protein: Mammalian Ste20-like protein kinase 4; MST4; RP6-213H19.1; Serine/threonine-protein kinase 26; STK26

## Application Notes

Optimal dilution of the MASK antibody should be determined by the researcher.

## Immunogen

Amino acids KLIEKFQKCSADESC were used as the immunogen for this MASK antibody.

## Storage

Aliquot and store the MASK antibody at -20oC.