


## MAST3 Antibody (R33627)

Catalog No.	Formulation	Size
R33627-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Predicted Reactivity	Cow, Dog
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	23031
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This MAST3 antibody is available for research use only.



Western blot testing of rat kidney lysate with MAST3 antibody at 0.3ug/ml. Predicted molecular weight: ~143kDa.

## Description

Additional name(s) for this target protein: Microtubule associated serine/threonine kinase 3

## Application Notes

Optimal dilution of the MAST3 antibody should be determined by the researcher.

## **Immunogen**

Amino acids DSFKKQEAVQEVS were used as the immunogen for this MAST3 antibody.

## **Storage**

Aliquot and store the MAST3 antibody at -20oC.