


Mboat4 Antibody (R33247)

| Catalog No. | Formulation | Size |
|--------------|---|--------|
| R33247-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

| | |
|-----------------------------|--|
| Availability | 1-3 business days |
| Species Reactivity | Human, Mouse |
| Predicted Reactivity | Cow, Pig |
| Format | Antigen affinity purified |
| Host | Goat |
| Clonality | Polyclonal (goat origin) |
| Isotype | Goat Ig |
| Purity | Antigen affinity |
| Gene ID | 619373 |
| Applications | Western Blot : 0.5-2ug/ml ELISA (peptide) LOD : 1:16000 |
| Limitations | This Mboat4 antibody is available for research use only. |



Western blot testing of mouse pancreas lysate with Mboat4 antibody at 0.5ug/ml. The expected ~50kDa band and the additional ~38kDa band are both blocked by the immunizing peptide.

Description

Additional name(s) for this target protein: Membrane bound O-acyltransferase domain containing 4

Application Notes

Optimal dilution of the Mboat4 antibody should be determined by the researcher.

Immunogen

Amino acids CHLGLHYTEYYLGEP were used as the immunogen for this Mboat4 antibody.

Storage

Aliquot and store the Mboat4 antibody at -20oC.