


MBP tag Antibody (F53952)

Catalog No.	Formulation	Size
F53952-0.2ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.2 ml
F53952-0.05ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.05 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	MBP tag
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	P0AEX9
Applications	Western Blot : 1:8000-16000
Limitations	This MBP tag antibody is available for research use only.



Western blot testing of E. coli lysate with MBP tag antibody at 1:16000. Predicted molecular weight: 43 kDa.



Western blot testing of MBP-tagged MART1 lysate with MBP tag antibody at 1:2000. Predicted molecular weight: 43 kDa.

Description

Involved in the high-affinity maltose membrane transport system MalEFGK. Initial receptor for the active transport of and chemotaxis toward maltooligosaccharides.

Application Notes

Titration of the MBP tag antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

MBP tag was used as the immunogen for this MBP tag antibody.

Storage

Aliquot the MBP tag antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.