


MC2R Antibody (R33892)

Catalog No.	Formulation	Size
R33892-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	4158
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This MC2R antibody is available for research use only.



Western blot testing of human skin lysate with MC2R antibody at 1ug/ml. Predicted molecular weight: ~33/43kDa (unmodified/glycosylated).

Description

Additional name(s) for this target protein: Melanocortin 2 receptor, ACTHR

Application Notes

Optimal dilution of the MC2R antibody should be determined by the researcher.

Immunogen

Amino acids ENINNTARNNSDC were used as the immunogen for this MC2R antibody.

Storage

Aliquot and store the MC2R antibody at -20oC.