


MC3R Antibody (R35931)

Catalog No.	Formulation	Size
R35931-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	4159
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This MC3R antibody is available for research use only.



Western blot testing of human brain lysate with MC3R antibody at 1ug/ml.
Predicted/observed molecular weight: ~40kDa.

Description

Additional name(s) for this target protein: Melanocortin 3 Receptor

Application Notes

Optimal dilution of the MC3R antibody should be determined by the researcher.

Immunogen

Amino acids SIQKTYLEGDFV were used as the immunogen for this MC3R antibody.

Storage

Aliquot and store the MC3R antibody at -20oC.