


MCTS1 Antibody (R33417)

Catalog No.	Formulation	Size
R33417-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Pig
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	28985
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This MCTS1 antibody is available for research use only.



Western blot testing of MOLT4 lysate with MCTS1 antibody at 0.3ug/ml. Predicted molecular weight: ~20kDa.

Description

Additional name(s) for this target protein: Malignant T cell amplified sequence 1

Application Notes

Optimal dilution of the MCTS1 antibody should be determined by the researcher.

Immunogen

Amino acids DEKENVSNCIQLKTS were used as the immunogen for this MCTS1 antibody.

Storage

Aliquot and store the MCTS1 antibody at -20oC.