


MCUR1 Antibody / Mitochondrial calcium uniporter regulator 1 (RQ4944)


Catalog No.	Formulation	Size
RQ4944	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

Bulk quote request


Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit IgG
Purity	Antigen affinity purified
Buffer	Lyophilized from 1X PBS with 2% Trehalose and 0.025% sodium azide
UniProt	Q96AQ8
Localization	Cytoplasmic
Applications	Western Blot : 0.5-1ug/ml Immunohistochemistry (FFPE) : 1-2ug/ml
Limitations	This MCUR1 antibody is available for research use only.




Western blot testing of human 1) HeLa, 2) MDA-MB-231, 3) HL-60, 4) MDA-MB-453, 5) A431, 6) Caco-2, 7) rat spleen, 8) mouse lung and 9) mouse ANA-1 lysate with MCUR1 antibody at 0.5ug/ml. Predicted molecular weight ~40 kDa.




IHC staining of FFPE human tonsil tissue with MCUR1 antibody at 1ug/ml. HIER: boil tissue sections in pH6, 10mM citrate buffer, for 10-20 min and allow to cool before testing.




IHC staining of FFPE human placenta tissue with MCUR1 antibody at 1ug/ml. HIER: boil tissue sections in pH6, 10mM citrate buffer, for 10-20 min and allow to cool before testing.



IHC staining of FFPE human lung cancer tissue with MCUR1 antibody at 1ug/ml. HIER: boil tissue sections in pH6, 10mM citrate buffer, for 10-20 min and allow to cool before testing.



IHC staining of FFPE human esophageal squamous cancer with MCUR1 antibody at 1ug/ml. HIER: boil tissue sections in pH6, 10mM citrate buffer, for 10-20 min and allow to cool before testing.



IHC staining of FFPE human ovarian cancer tissue with MCUR1 antibody at 1ug/ml. HIER: boil tissue sections in pH6, 10mM citrate buffer, for 10-20 min and allow to cool before testing.

Description

MCUR1 is an inner mitochondrial membrane protein that has been shown to function as a component of the mitochondrial $\text{Ca}(2+)$ uniporter, in the assembly of mitochondrial respiratory complex IV, and in mitochondrial permeability transition (MPT). The MCUR1 gene is mapped to chromosome 6p23 based on an alignment of the MCUR1 sequence with the genomic sequence.

Application Notes

Optimal dilution of the MCUR1 antibody should be determined by the researcher.

Immunogen

Amino acids ATQQAELIVSALVKILEANMDIVYKDMVTKMQQE from the human protein were used as the immunogen for the MCUR1 antibody.

Storage

After reconstitution, the MCUR1 antibody can be stored for up to one month at 4°C. For long-term, aliquot and store at

-20oC. Avoid repeated freezing and thawing.