


MEIS1 Antibody [clone ACDO-13] (RQ4694)

Catalog No.	Formulation	Size
RQ4694	Antibody in PBS with 0.02% sodium azide, 50% glycerol and 0.4-0.5mg/ml BSA	100 ul

[Bulk quote request](#)

Availability	1-2 weeks
Species Reactivity	Human
Format	Purified
Host	Rabbit
Clonality	Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	ACDO-13
Purity	Affinity purified
UniProt	O00470
Applications	Western Blot : 1:1000-1:5000
Limitations	This MEIS1 antibody is available for research use only.



Western blot testing of human K562 cell lysate with MEIS1 antibody. Predicted molecular weight: ~43 kDa.

Description

Homeobox protein Meis1 is a MEIS protein belonging to the TALE (Three Amino Acid Loop Extension) homeobox containing transcription factor family. MEIS1 has been associated with leukemogenesis and neuroblastoma.

Application Notes

Optimal dilution of the MEIS1 antibody should be determined by the researcher.

Immunogen

A synthetic peptide specific to human MEIS1 was used as the immunogen for the MEIS1 antibody.

Storage

Store the MEIS1 antibody at -20oC.