


Membrane Palmitoylated Protein 1 Antibody / MPP1 [clone 28M87] (RQ8529)

Catalog No.	Formulation	Size
RQ8529	Antibody in PBS with 0.02% sodium azide, 50% glycerol and 0.4-0.5mg/ml BSA	100 ul

Recombinant **RABBIT MONOCLONAL**

[Bulk quote request](#)

Availability	1-3 days
Species Reactivity	Human
Format	Purified
Host	Rabbit
Clonality	Recombinant Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	28M87
Purity	Affinity chromatography
UniProt	Q00013
Applications	Western Blot : 1:500-1:2000
Limitations	This Membrane Palmitoylated Protein 1 antibody is available for research use only.



Western blot testing of human Jurkat cell lysate with Membrane Palmitoylated Protein 1 antibody. Expected molecular weight 52-55 kDa.

Description

MPP1 / Membrane Palmitoylated Protein 1, also called 55 kDa erythrocyte membrane protein, is an essential regulator of neutrophil polarity. [UniProt]

Application Notes

Optimal dilution of the Membrane Palmitoylated Protein 1 antibody should be determined by the researcher.

Immunogen

A synthetic peptide from human Membrane Palmitoylated Protein 1 protein was used as the immunogen for the Membrane Palmitoylated Protein 1 antibody.

Storage

Store the Membrane Palmitoylated Protein 1 antibody at -20oC.