


MEX3C Antibody (R36247)

Catalog No.	Formulation	Size
R36247-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Cow, Dog, Pig, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	51320
Applications	Western Blot : 0.2-0.5ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This MEX3C antibody is available for research use only.



Western blot testing of mouse testis lysate with MEX3C antibody at 0.2ug/ml. The ~63/36kDa bands are significantly reduced.

Description

Additional name(s) for this target protein: Mex-3 homolog C

Application Notes

Optimal dilution of the MEX3C antibody should be determined by the researcher.

Immunogen

Amino acids TPRLSPTFPESIEH were used as the immunogen for this MEX3C antibody.

Storage

Aliquot and store the MEX3C antibody at -20oC.