


MICS1 Antibody (R35160)


Catalog No.	Formulation	Size
R35160-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat, Dog, Pig, Cow
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	27069
Applications	Western Blot : 0.01-0.1ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This MICS1 antibody is available for research use only.



Western blot testing of HeLa cell lysate with MICS1 antibody at 0.3ug/ml. Predicted molecular weight: ~38 kDa.



Western blot testing of human cerebellum (A), duodenum (B) and heart (C) lysates with MICS1 antibody at 0.01ug/ml. Predicted molecular weight: ~38 kDa.

Description

Additional name(s) for this target protein: Mitochondrial morphology and cristae structure 1, Growth hormone inducible transmembrane protein, GHITM

Application Notes

Optimal dilution of the MICS1 antibody should be determined by the researcher.

Immunogen

Amino acids PSREYATKTRIGIR were used as the immunogen for this MICS1 antibody.

Storage

Aliquot and store the MICS1 antibody at -20oC.