


MKRN1 Antibody (R35818)


Catalog No.	Formulation	Size
R35818-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)


Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	23608
Localization	Nuclear, cytoplasmic
Applications	Western Blot : 0.3-1ug/ml IF/ICC : 10ug/ml Flow Cytometry : 10ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This MKRN1 antibody is available for research use only.




IF/ICC staining of fixed and permeabilized human U-2 OS cells with MKRN1 antibody (green) at 10ug/ml and DAPI nuclear stain (blue).



Western blot testing of human testis tissue lysate with MKRN1 antibody. Predicted molecular weight ~53 kDa.



Western blot testing of human ThP-1 cell lysate with MKRN1 antibody. Predicted molecular weight ~53 kDa.



FACS testing of fixed and permeabilized human ThP-1 cells with MKRN1 antibody (blue) at 10ug/ml and naive goat Ig (black).

Description

Additional name(s) for this target protein: Makorin ring finger protein 1

Application Notes

Optimal dilution of the MKRN1 antibody should be determined by the researcher.

Immunogen

Amino acids RYEHSKPLKQEEAT were used as the immunogen for this MKRN1 antibody.

Storage

Aliquot and store the MKRN1 antibody at -20oC.