


## MKRN1 Antibody (R36342)

Catalog No.	Formulation	Size
R36342-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	23608
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This MKRN1 antibody is available for research use only.



Western blot: transfected and FLAG tagged HEK293 lysate probed with Lane A: MKRN1 antibody [1ug/ml], Lane C: anti-FLAG [1/1000] and mock-transfected HEK293 probed with MKRN1 antibody [1ug/ml] in Lane B. Observed molecular weight: 50~60kDa.

## Description

Additional name(s) for this target protein: Makorin

## Application Notes

Optimal dilution of the MKRN1 antibody should be determined by the researcher.

1. This MKRN1 antibody may cross react with MKRN4 (~53 kDa).

## Immunogen

Amino acids VTKEESEKEQTAVET were used as the immunogen for this MKRN1 antibody.

## Storage

Aliquot and store the MKRN1 antibody at -20oC.