


MMP14 Antibody (R33985)

Catalog No.	Formulation	Size
R33985-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Pig, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	4323
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This MMP14 antibody is available for research use only.



Western blot testing of mouse NIH3T3 [A], A549 [B], and HeLa [C] lysates with MMP14 antibody at 1ug/ml. Predicted molecular weight: ~66/54kDa.

Description

Additional name(s) for this target protein: Matrix metalloproteinase 14

Application Notes

Optimal dilution of the MMP14 antibody should be determined by the researcher.

Immunogen

Amino acids ESATPLRFREVP were used as the immunogen for this MMP14 antibody.

Storage

Aliquot and store the MMP14 antibody at -20oC.