


## MMP3 Antibody (R33577)


Catalog No.	Formulation	Size
R33577-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

<b>Availability</b>	1-3 business days
<b>Species Reactivity</b>	Human, Mouse
<b>Format</b>	Antigen affinity purified
<b>Host</b>	Goat
<b>Clonality</b>	Polyclonal (goat origin)
<b>Isotype</b>	Goat Ig
<b>Purity</b>	Antigen affinity
<b>Gene ID</b>	4314
<b>Applications</b>	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:128000
<b>Limitations</b>	This MMP3 antibody is available for research use only.



Western blot testing of human prostate lysate with MMP3 antibody at 0.2ug/ml. Molecular weight: the expected ~54 kDa band and the additional ~24kDa band are both blocked by the immunizing peptide.



Western blot testing of mouse eye lysate with MMP3 antibody at 0.3ug/ml. Predicted molecular weight ~54 kDa.

## Description

Additional name(s) for this target protein: Matrix metalloproteinase 3

## Application Notes

Optimal dilution of the MMP3 antibody should be determined by the researcher.

## Immunogen

Amino acids DPNAKKVTHTLKSNS were used as the immunogen for this MMP3 antibody.

## Storage

Aliquot and store the MMP3 antibody at -20oC.