


MOGAT2 Antibody (R34417)

Catalog No.	Formulation	Size
R34417-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Pig
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	80168
Applications	Western Blot : 0.05-0.2ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This MOGAT2 antibody is available for research use only.



Western blot testing of human duodenum lysate with MOGAT2 antibody at 0.05ug/ml.
Predicted/observed molecular weight: ~37kDa.

Description

Additional name(s) for this target protein: Monoacylglycerol O-acyltransferase 2

Application Notes

Optimal dilution of the MOGAT2 antibody should be determined by the researcher.

Immunogen

Amino acids KESAAHILNRK were used as the immunogen for this MOGAT2 antibody.

Storage

Aliquot and store the MOGAT2 antibody at -20oC.