


MORF4L2 Antibody (R36074)


Catalog No.	Formulation	Size
R36074-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	9643
Applications	Western Blot : 0.2-1ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This MORF4L2 antibody is available for research use only.



Western blot testing of A431 lysate with MORF4L2 antibody at 0.2ug/ml. Predicted molecular weight: ~32kDa.



Western blot of HEK293 lysate overexpressing MORF4L2 probed with MORF4L2 antibody (mock transfection in lane 1). Predicted molecular weight: ~32kDa.

Description

Additional name(s) for this target protein: MORF4L2, MRGX, MORF-related gene X

Application Notes

Optimal dilution of the MORF4L2 antibody should be determined by the researcher.

Immunogen

Amino acids SSRKQGSQPRGQQS were used as the immunogen for this MORF4L2 antibody.

Storage

Aliquot and store the MORF4L2 antibody at -20oC.