


## MORF4L2 Antibody (R36074)


Catalog No.	Formulation	Size
R36074-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	9643
Applications	Western Blot : 0.2-1ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This MORF4L2 antibody is available for research use only.



Western blot testing of A431 lysate with MORF4L2 antibody at 0.2ug/ml. Predicted molecular weight: ~32kDa.



Western blot of HEK293 lysate overexpressing MORF4L2 probed with MORF4L2 antibody (mock transfection in lane 1). Predicted molecular weight: ~32kDa.

## Description

Additional name(s) for this target protein: MORF4L2, MRGX, MORF-related gene X

## Application Notes

Optimal dilution of the MORF4L2 antibody should be determined by the researcher.

## Immunogen

Amino acids SSRKQGSQPRGQQS were used as the immunogen for this MORF4L2 antibody.

## Storage

Aliquot and store the MORF4L2 antibody at -20oC.