

Mouse Interferon gamma Antibody / IFNG (RQ4622)

Catalog No.	Formulation	Size
RQ4622	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Mouse
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit IgG
Purity	Antigen affinity purified
Buffer	Lyophilized from 1X PBS with 2% Trehalose and 0.025% sodium azide
UniProt	P01580
Applications	ELISA (Capture) : 2-8ug/ml (mouse recombinant protein)
Limitations	This mouse Interferon gamma antibody is available for research use only.



Description

Interferon-gamma (IFN-gamma) is an inflammatory cytokine that has been implicated in the development of fibrosis in inflamed tissues. The production of IFN-gamma, which is under genetic control, can influence the development of fibrosis in lung allografts. IFN-gamma is also produced by natural killer (NK) cells and most prominently by CD8 cytotoxic T cells, and is vital for the control of microbial pathogens. Interferon gamma is believed to be crucial for host defence against many infections. Genetically determined variability in IFN-gamma and expression might be important for the development of tuberculosis. IFN-gamma activates human macrophage oxidative metabolism and antimicrobial activity. In addition to having antiviral activity, IFN-gamma has important immunoregulatory functions. IFN-gamma plays an important role in the

control of neointima proliferation.

Application Notes

Optimal dilution of the mouse Interferon gamma antibody should be determined by the researcher.

Immunogen

Amino acids M1-C155 from the mouse protein were used as the immunogen for the mouse Interferon gamma antibody.

Storage

After reconstitution, the mouse Interferon gamma antibody can be stored for up to one month at 4oC. For long-term, aliquot and store at -20oC. Avoid repeated freezing and thawing.