


## MPG Antibody (R36105)

Catalog No.	Formulation	Size
R36105-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	4350
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:1000
Limitations	This MPG antibody is available for research use only.



Western blot testing of HEK293 lysate with MPG antibody at 1ug/ml. Predicted molecular weight: ~33kDa.

## Description

Additional name(s) for this target protein: N-methylpurine-DNA glycosylase (MPG), N-methylpurine-DNA glycosylase (MDG)

## Application Notes

Optimal dilution of the MPG antibody should be determined by the researcher.

## **Immunogen**

Amino acids SVVDRVAEQDTQA were used as the immunogen for this MPG antibody.

## **Storage**

Aliquot and store the MPG antibody at -20oC.