


MPG Antibody (R36105)

Catalog No.	Formulation	Size
R36105-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	4350
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:1000
Limitations	This MPG antibody is available for research use only.



Western blot testing of HEK293 lysate with MPG antibody at 1ug/ml. Predicted molecular weight: ~33kDa.

Description

Additional name(s) for this target protein: N-methylpurine-DNA glycosylase (MPG), N-methylpurine-DNA glycosylase (MDG)

Application Notes

Optimal dilution of the MPG antibody should be determined by the researcher.

Immunogen

Amino acids SVVDRVAEQDTQA were used as the immunogen for this MPG antibody.

Storage

Aliquot and store the MPG antibody at -20oC.