


MPP5 Antibody (R34832)

Catalog No.	Formulation	Size
R34832-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Pig
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	64398
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This MPP5 antibody is available for research use only.



Western blot testing of human brain lysate with MPP5 antibody at 1ug/ml. Predicted molecular weight: ~77kDa.

Description

Additional name(s) for this target protein: Membrane protein palmitoylated 5, PALS, PALS1; MPP5

Application Notes

Optimal dilution of the MPP5 antibody should be determined by the researcher.

Immunogen

Amino acids NKASPPFPLISNAQ were used as the immunogen for this MPP5 antibody.

Storage

Aliquot and store the MPP5 antibody at -20oC.