


MRGBP Antibody (R33253)

Catalog No.	Formulation	Size
R33253-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	55257
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This MRGBP antibody is available for research use only.



Western blot testing of human testis lysate with MRGBP antibody at 0.3ug/ml. Predicted molecular weight: ~22kDa.

Description

Additional name(s) for this target protein: MRG(MORF4-related gene)-binding protein

Application Notes

Optimal dilution of the MRGBP antibody should be determined by the researcher.

Immunogen

Amino acids CIRDKFSQNIGRQVP were used as the immunogen for this MRGBP antibody.

Storage

Aliquot and store the MRGBP antibody at -20oC.