


## MRP8 Antibody / ABCC11 (R34694)


Catalog No.	Formulation	Size
R34694-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
UniProt	Q96J66
Gene ID	85320
Applications	Western Blot : 0.5-1.5ug/ml IHC (FFPE) : 3.75ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This MRP8 antibody is available for research use only.



Western blot testing of human breast lysate with MRP8 antibody at 0.5ug/ml. Predicted molecular weight: ~154/150 kDa (isoforms 1/2), observed here at ~125 kDa.



IHC testing of FFPE human breast tissue with MRP8 antibody at 3.75ug/ml. Required HIER: steamed antigen retrieval with pH6 citrate buffer; AP-staining.

## Description

Additional name(s) for this target protein: Multi-resistance protein 8, ATP-binding cassette, sub-family C (CFTR/MRP), member 11; ABCC11

## Application Notes

Optimal dilution of the MRP8 antibody should be determined by the researcher.

## Immunogen

Amino acids KVVEFDRPEVLRK were used as the immunogen for this MRP8 antibody.

## Storage

Aliquot and store the MRP8 antibody at -20oC.