


MSC Antibody (R35056)

Catalog No.	Formulation	Size
R35056-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	9242
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This MSC antibody is available for research use only.



Western blot testing of human bone marrow lysate with MSC antibody at 1ug/ml.
Predicted molecular weight: ~24kDa.

Description

Additional name(s) for this target protein: Musculin, Activated B-cell factor-1, ABF1

Application Notes

Optimal dilution of the MSC antibody should be determined by the researcher.

Immunogen

Amino acids PLKSDWNSLDGIIR were used as the immunogen for this MSC antibody.

Storage

Aliquot and store the MSC antibody at -20oC.