

MSH2 Antibody / MutS homolog 2 [clone MSVA-902M] (V6099)

Catalog No.	Formulation	Size
V6099-100UG	Antibody in 1X PBS with 0.05% BSA, 0.05% sodium azide	100 ug
V6099-20UG	Antibody in 1X PBS with 0.05% BSA, 0.05% sodium azide	20 ug

Recombinant **MOUSE MONOCLONAL**

[Bulk quote request](#)

Species Reactivity	Human
Format	Purified
Host	Mouse
Clonality	Recombinant Mouse Monoclonal
Isotype	Mouse IgG2a, kappa
Clone Name	MSVA-902M
UniProt	IDP43246
Localization	Chromosome, Nucleus
Applications	Immunohistochemistry (FFPE) : 1-2ug/ml
Limitations	This MSH2/MutS homolog 2 antibody is available for research use only.

Application Notes

1. Optimal dilution of the MSH2/MutS homolog 2 antibody should be determined by the researcher.
2. This MSH2/MutS homolog 2 antibody is recombinantly produced by expression in CHO cells.

Immunogen

Recombinant full-length human MSH2 protein was used as the immunogen for the MSH2/MutS homolog 2 antibody.

Storage

MSH2/MutS homolog 2 antibody with sodium azide - store at 2 to 8°C; antibody without sodium azide - store at -20 to -80°C.

