


MSH3 Antibody [clone RM405] (R20421)

Catalog No.	Formulation	Size
R20421-0.1ML	Antibody in PBS with 50% glycerol, 1% BSA and 0.09% sodium azide	100 ul


Recombinant **RABBIT MONOCLONAL**

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Host	Rabbit
Clonality	Recombinant Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	RM405
Purity	Protein A purified from animal origin-free supernatant
UniProt	P20585
Localization	Nuclear
Applications	Immunohistochemistry (FFPE) : 1:100-1:200 Western Blot : 1:400-1:800
Limitations	This recombinant MSH3 antibody is available for research use only.



Western blot testing of human HEK293 cell lysate with recombinant MSH3 antibody at 1:800 dilution. Predicted molecular weight ~127 kDa.



IHC staining of FFPE human colon cancer tissue with recombinant MSH3 antibody at 1:200.

Description

This antibody reacts to human to MSH3 / DNA mismatch repair protein Msh3.

Application Notes

The stated application concentrations are suggested starting points. Titration of the recombinant MSH3 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A peptide corresponding to residues near the N-terminus of human MSH3 was used as the immunogen for the recombinant MSH3 antibody.

Storage

Store the recombinant MSH3 antibody at -20oC.