


MSH5 Antibody (R34074)

Catalog No.	Formulation	Size
R34074-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	4439
Applications	Immunohistochemistry (FFPE) : 2.5-5ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This MSH5 antibody is available for research use only.



IHC testing of FFPE human testis with MSH5 antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

Description

Additional name(s) for this target protein: MutS homolog 5

Application Notes

Optimal dilution of the MSH5 antibody should be determined by the researcher.

Immunogen

Amino acids GDLHCEIRDQET were used as the immunogen for this MSH5 antibody.

Storage

Aliquot and store the MSH5 antibody at -20oC.