


MST2 Antibody (R33439)

Catalog No.	Formulation	Size
R33439-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Cow, Dog, Pig, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
UniProt	Q13188
Gene ID	6788
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This MST2 antibody is available for research use only.



Western blot testing of NIH3T3 lysate with MST2 antibody at 1ug/ml. Predicted molecular weight: ~59/56kDa (isoforms 1/2).

Description

Additional name(s) for this target protein: Serine/threonine kinase 3; STK3

Application Notes

Optimal dilution of the MST2 antibody should be determined by the researcher.

Immunogen

Amino acids DFKNKSHENCNQN were used as the immunogen for this MST2 antibody.

Storage

Aliquot and store the MST2 antibody at -20oC.