


MT-ATP6 Antibody (R35171)

Catalog No.	Formulation	Size
R35171-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	4508
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This MT-ATP6 antibody is available for research use only.



Western blot testing of human heart lysate with MT-ATP6 antibody at 0.1ug/ml.
Predicted molecular weight: ~25kDa.

Description

Additional name(s) for this target protein: F-ATPase protein 6, ATPase6; ATP6

Application Notes

Optimal dilution of the MT-ATP6 antibody should be determined by the researcher.

Immunogen

Amino acids PTSKYLINNRLITTQ were used as the immunogen for this MT-ATP6 antibody.

Storage

Aliquot and store the MT-ATP6 antibody at -20oC.