


## MTHFS Antibody (R34607)

Catalog No.	Formulation	Size
R34607-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	10588
Applications	Western Blot : 0.5-1.5ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This MTHFS antibody is available for research use only.



Western blot testing of human liver lysate with MTHFS antibody at 0.5ug/ml with [B] and without [A] blocking/immunizing peptide. Predicted molecular weight: ~23kDa.

## Description

Additional name(s) for this target protein: 5,10-methenyltetrahydrofolate synthetase

## Application Notes

Optimal dilution of the MTHFS antibody should be determined by the researcher.

## Immunogen

Amino acids KRCLQHGEVKPYT were used as the immunogen for this MTHFS antibody.

## Storage

Aliquot and store the MTHFS antibody at -20oC.