


## MUC16 Antibody / Cancer antigen 125 [clone CA125/9923R] (V6041)

Catalog No.	Formulation	Size
V6041-100UG	0.2 mg/ml in 1X PBS with 0.05% BSA, 0.05% sodium azide	100 ug
V6041-20UG	0.2 mg/ml in 1X PBS with 0.05% BSA, 0.05% sodium azide	20 ug
V6041SAF-100UG	1 mg/ml in 1X PBS; BSA free, sodium azide free	100 ug


Recombinant **RABBIT MONOCLONAL**

[Bulk quote request](#)


<b>Species Reactivity</b>	Human
<b>Format</b>	Purified
<b>Host</b>	Rabbit
<b>Clonality</b>	Recombinant Rabbit Monoclonal
<b>Isotype</b>	Rabbit IgG, kappa
<b>Clone Name</b>	CA125/9923R
<b>UniProt</b>	Q8WXI7
<b>Localization</b>	Cell membrane, Extracellular space, Secreted
<b>Applications</b>	ELISA : Immunohistochemistry (FFPE) : 1-2ug/ml
<b>Limitations</b>	This MUC16/Cancer antigen 125 antibody is available for research use only.



Formalin-fixed, paraffin-embedded human colon carcinoma stained with MUC16 Recombinant Rabbit Monoclonal Antibody (CA125/9923R). Inset: PBS instead of primary antibody; secondary only negative control.



SDS-PAGE Analysis of Purified MUC16 Recombinant Rabbit Monoclonal Antibody (CA125/9923R). Confirmation of Purity and Integrity of Antibody.



Formalin-fixed, paraffin-embedded human ovarian carcinoma stained with MUC16 Recombinant Rabbit Monoclonal Antibody (CA125/9923R). Inset: PBS instead of primary antibody; secondary only negative control.

## Application Notes

Optimal dilution of the MUC16/Cancer antigen 125 antibody should be determined by the researcher.

## Immunogen

Recombinant human MUC16 protein was used as the immunogen for the MUC16/Cancer antigen 125 antibody.

## Storage

MUC16/Cancer antigen 125 antibody with sodium azide - store at 2 to 8°C; antibody without sodium azide - store at -20 to -80°C.