



## MUC3B Antibody (F41684)

Catalog No.	Formulation	Size
F41684-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F41684-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

<b>Availability</b>	1-3 business days
<b>Species Reactivity</b>	Human
<b>Format</b>	Antigen affinity purified
<b>Host</b>	Rabbit
<b>Clonality</b>	Polyclonal (rabbit origin)
<b>Isotype</b>	Rabbit Ig
<b>Purity</b>	Antigen affinity
<b>UniProt</b>	Q9H195
<b>Applications</b>	Western Blot : 1:1000 IHC (Paraffin) : 1:10-1:50
<b>Limitations</b>	This MUC3B antibody is available for research use only.



## Description

Major glycoprotein component of a variety of mucus gels. Thought to provide a protective, lubricating barrier against particles and infectious agents at mucosal surfaces (By similarity).

## Application Notes

Titration of the MUC3B antibody may be required due to differences in protocols and secondary/substrate sensitivity.

## Immunogen

A portion of amino acids 866-895 from the human protein was used as the immunogen for this MUC3B antibody.

## Storage

Aliquot the MUC3B antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.