

MUC6 Antibody / Mucin 6 [clone MSVA-806R] (V5954)

Catalog No.	Formulation	Size
V5954-100UG	0.2 mg/ml in 1X PBS with 0.05% BSA, 0.05% sodium azide	100 ug
V5954-20UG	0.2 mg/ml in 1X PBS with 0.05% BSA, 0.05% sodium azide	20 ug

Recombinant **RABBIT MONOCLONAL**

[Bulk quote request](#)

Species Reactivity	Human
Format	Purified
Host	Rabbit
Clonality	Recombinant Rabbit Monoclonal
Isotype	Rabbit IgG, kappa
Clone Name	MSVA-806R
UniProt	Q6W4X9
Localization	Secreted
Applications	Immunohistochemistry (FFPE) : 1-2ug/ml
Limitations	This MUC6/Mucin 6 antibody is available for research use only.

Application Notes

1. Optimal dilution of the MUC6/Mucin 6 antibody should be determined by the researcher.
2. This MUC6/Mucin 6 antibody is recombinantly produced by expression in human HEK293 cells.

Immunogen

Recombinant human MUC6 protein was used as the immunogen for the MUC6/Mucin 6 antibody.

Storage

MUC6/Mucin 6 antibody with sodium azide - store at 2 to 8°C; antibody without sodium azide - store at -20 to -80°C.

