


MUNC18-1 Antibody / STXBP1 (R36062)


Catalog No.	Formulation	Size
R36062-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	6812
Applications	Western Blot : 0.1-0.3ug/ml Immunohistochemistry (FFPE) : 5-10ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This MUNC18-1 antibody is available for research use only.



Western blot testing of human brain lysate with MUNC18-1 antibody at 0.1ug/ml.
Predicted molecular weight: ~69kDa.



IHC staining of FFPE human cortex with MUNC18-1 antibody at 5ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

Description

Additional name(s) for this target protein: Syntaxin binding protein 1, STXBP1

Application Notes

Optimal dilution of the MUNC18-1 antibody should be determined by the researcher.

Immunogen

Amino acids SRVSFEDQAPTME were used as the immunogen for this MUNC18-1 antibody.

Storage

Aliquot and store the MUNC18-1 antibody at -20oC.