


## MUNC18-3 Antibody (R35945)

| Catalog No.  | Formulation   | Size   |
|--------------|---|--------|
| R35945-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

|                      |   |
|----------------------|---|
| Availability         | 1-3 business days   |
| Species Reactivity   | Human, Mouse, Rat   |
| Predicted Reactivity | Dog, Pig  |
| Format               | Antigen affinity purified                                   |
| Host                 | Goat  |
| Clonality            | Polyclonal (goat origin)                                    |
| Isotype              | Goat Ig   |
| Purity               | Antigen affinity  |
| Gene ID              | 6814  |
| Applications         | Western Blot : 0.3-1ug/ml<br>ELISA (peptide) LOD : 1:128000 |
| Limitations          | This MUNC18-3 antibody is available for research use only.  |



Western blot testing of mouse [A] and rat [B] spleen lysates with MUNC18-3 antibody at 1ug/ml. Predicted molecular weight: ~68kDa.

### Description

Additional name(s) for this target protein: STXBP3, syntaxin binding protein 3, MUNC18C

### Application Notes

Optimal dilution of the MUNC18-3 antibody should be determined by the researcher.

## Immunogen

Amino acids SKSIRRCKEINIS were used as the immunogen for this MUNC18-3 antibody.

## Storage

Aliquot and store the MUNC18-3 antibody at -20oC.