


MUNC18C Antibody (R35732)

Catalog No.	Formulation	Size
R35732-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Pig
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	6814
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This MUNC18C antibody is available for research use only.



Western blot testing of human platelet lysate with MUNC18C antibody at 0.1ug/ml.
Predicted molecular weight: ~68kDa.

Description

Additional name(s) for this target protein: STXBP3, syntaxin binding protein 3, MUNC18-3

Application Notes

Optimal dilution of the MUNC18C antibody should be determined by the researcher.

Immunogen

Amino acids RQKPRANYLEDNRKN were used as the immunogen for this MUNC18C antibody.

Storage

Aliquot and store the MUNC18C antibody at -20oC.