


MYB Antibody (R34267)

Catalog No.	Formulation	Size
R34267-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Pig, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	4602
Applications	Western Blot : 0.5-1.5ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This MYB antibody is available for research use only.



Western blot testing of Jurkat lysate with MYB antibody at 0.5ug/ml. Predicted molecular weight 37~85 kDa (isoforms 1-12), observed here at ~85kDa..

Description

Additional name(s) for this target protein: C-Myb, v-Myb

Application Notes

Optimal dilution of the MYB antibody should be determined by the researcher.

1. The amino acid sequence used as immunogen is common to all 12 isoforms.

Immunogen

Amino acids HRDQTVKTQKENT were used as the immunogen for this MYB antibody. (1)

Storage

Aliquot and store the MYB antibody at -20oC.