


## Myelin Protein Zero Antibody / MPZ (R33475)

Catalog No.	Formulation	Size
R33475-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Predicted Reactivity	Dog
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	4359
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This Myelin Protein Zero antibody is available for research use only.



Western blot testing of rat spinal cord lysate with Myelin Protein Zero antibody at 0.1ug/ml. Predicted molecular weight: ~28kDa.

### Description

Additional name(s) for this target protein: MPZ

### Application Notes

Optimal dilution of the Myelin Protein Zero antibody should be determined by the researcher.

## **Immunogen**

Amino acids DHSRSTKAVSEK were used as the immunogen for this Myelin Protein Zero antibody.

## **Storage**

Aliquot and store the Myelin Protein Zero antibody at -20oC.